


AUTOMATIC STAINER AUTOCOLOR TOUCH


Code: 16-1600/T



USE PRECAUTIONS

Before using the instrument, read carefully the instructions and warnings contained in this manual and keep it for further reference. They supply important indications regarding the functions and safety for installing, using and maintaining the instrument.

Bio-Optica Milano S.p.A. cannot be held responsible for any damage caused by improper or incorrect use and by the non-observance of any of the prescription provided in this manual and by the safety regulations in force.

1. After unpacking, make sure that the instrument is complete and not damaged by transport.
 2. Before connecting the instrument to the power supply make sure that its rating corresponds to that of the power supply.
 3. This instrument must only be used for the purpose for which it was designed, that is, as automatic stainer for laboratory use. Any other use is to be considered improper and therefore hazardous.
 4. The instrument must only be used by authorized and professionally qualified technician.
 5. The electrical safety of this instrument can be guaranteed only if it is correctly connected to an efficient earth circuit as indicated by current electrical safety regulations. It is necessary to check this fundamental safety prerequisite, and if in doubt, ask to check the circuit. The instrument is provided with a power supply cable having 2 wires + ground tap that have to be connected to the power supply socket.
 6. Do not remove the chassis or parts of it during operation. Switch off the instrument and disconnect the power supply cable before opening it. This operation must to be effected only by authorized and professionally qualified technician.
 7. To eliminate instrument malfunctioning risks, do not work near strong magnetic fields and do not use transmitters such as cellular phones near the instrument. In case of serious malfunctioning switch off the instrument and contact the Technical Assistance Service.
 8. All waste material, both infectious and radioactive, deriving from the appliance working cycle must be disposed in compliance with the regulation in force.
-  This appliance is marked from this symbol, in compliance with EU directive 2002/96/CE regarding electric and electronic appliances waste. This mean that the instrument, at the end of its useful life, must be collected separately from other refuse. The user must deliver it to the special differentiated refuse collection centres, that are predisposed by the public authority.
9. The contents of this manual is subject to change without further notice.
 10. Please find enclosed the declaration of conformity.

11. Graphic symbols indicated on the label (positioned near the instrument's power supply socket):

Symbol for CATALOGUE NUMBER:



Symbol for SERIAL NUMBER:



Symbol for ALTERNATING CURRENT:



Symbol for FUSE:



Symbol for CONSULT THE INSTRUCTIONS:



Symbol for EC MARK:



Symbol for IN VITRO DIAGNOSTIC-MEDICAL DEVICE:



Symbol for DISPOSAL OF ELECTRIC AND ELECTRONIC EQUIPMENT:



Symbol for DATE OF MANUFACTURE:



Symbol for MANUFACTURER:



INDEX

1) Technical features.	Page 5
2) Installation.	Page 7
3) Water connections.	Page 7
4) Control panel.	Page 10
5) Instrument's start.	Page 11
6) Date and hour setup.	Page 13
7) Dishes map and relative contents.	Page 14
8) Parameters setup and programming.	Page 14
9) Alarms.	Page 23
10) Maintenance.	Page 24
11) Filters.	Page 25
12) Main safety elements.	Page 25
13) Instructions for the rotation of the dish holder disk.	Page 26
14) Replacing the fuses.	Page 27
15) Cleaning and maintenance.	Page 27
16) Main components.	Page 28
17) Problem solving.	Page 30
18) Electric wiring diagram.	Page 31
19) Declaration of conformity.	Page 35

1) TECHNICAL FEATURES

Carousel automatic stainer for colouring histological and cytological slides.

Dimensional features

Dimensions (WxDxH): 600x700x600 mm.

Weight: 75 Kg.

Electrical connections

Power supply: 230V~ 50/60Hz.

Power: 120 Watt.

Other connections

Water connections: Provided pipes for supplying the running water dishes:

N. 1 running water inlet pipe, internal Ø 10 mm, provided with rapid connection with interception.

N. 1 running water outlet pipe, internal Ø 18 mm.

N. 2 neoprene pipes - internal Ø 10 mm, one metre long - provided with rapid connection with interception (for the connection between the dishes and the instrument).

Water pressure: From 2,5 to 4 bar.

Structural features

- Painted steel chassis resistant to scratches and solvents.
On the right side are placed: N. 1 green general switch, N. 1 power supply socket with relative fuses, N. 2 connectors for pipes connecting running water dishes/instrument, N. 1 warning light for line voltage existence.
On the left side are placed: N. 1 connector connecting running water inlet pipe, N. 1 connector connecting running water outlet pipe, N. 1 running water pressure regulator, N. 1 filter seat cover provided with handle.
- External cover for containing fumes.
- Basket holder cover, positioned under the external cover, provided with 24 guides for coupling one or more baskets and 24 stainless steel covers for covering the dishes.
- Central basin for carrying waste water from wash dishes to the drain placed on the left side of the instrument.
- N. 4 adjustable rubber feet to get a perfect horizontal position.
- Aspiration/filtration system for the treatment of fumes. The instrument is provided with N. 1 active charcoal cartridge filter specific for alcohol and xilol. Practical filter's replacement from the left side of the instrument.
- Revolving dish holder disk to facilitate the access to the dishes placed on the back side.
- N. 2 stainless steel slide holder baskets (capacity 65 standard slides 26x76 mm).

- N. 21 reagents dishes (capacity 750 ml) in white thermoplastic material.
 - N. 1 dish (capacity 750 ml) in white thermoplastic material for baskets load/discharge in position 0.
 - N. 2 running water dishes (capacity 750 ml) in transparent thermoplastic material.
 - N. 2 independent electric valves which activate the water flow when there is the passage in running water.
 - Pressure regulator to regulate the water flow.
 - N. 2 lithium buffer batteries for emergency power supply in case of black-out. Instrument automatic restart when returning the current.
-

Working features

- General ignition switch.
- Warning light of line voltage.
- Emergency switch for stopping all the machine functions at any time.
- Visual and acoustic alarm in case of: cycle end, filter's saturation, absence of current.
- Visual and acoustic alarm in case of table out of position.
- Visual and acoustic alarm for water discharge (malfunction of electric valve, obstruction discharge pipes, etc.).
- Page of visualization and description of alarms accessible also during the execution of the programme.
- Possibility of compiling the map of the reagents loaded in machine specifying: name of reagent, date of insertion, duration of reagent.
- Possibility of visualization of programme steps during the execution of the cycle.
- Possibility of saving up to 49 different staining protocols, with a maximum number of 30 passages each.
- Possibility of programming for each station: time of immersion and dripping from 00' : 01" up to 50' : 59".
- Possibility of setting up in machine parameters time of pause between the stations from 00' : 01" up to 50' : 59". Pause time is set up by default at 00' : 01".
- Possibility of changing, during the execution of the programme, the speed of oscillation, dripping and speed of rotation.
- Possibility of programming the filter's replacement.
- Possibility of stopping and restarting a programme (in case of possible controls, filling up of reagents, ecc.).
- Double reset of safety in case of accidental stop.
- Possibility of visualizing the map of the stations during the cycle.
- Entire management of the instrument through use of single touch screen panel.

2) INSTALLATION

During the stainer moving, proceed cautiously in order to avoid possible damages.

Position the instrument on a level and stable working bench, free of vibrations. Make sure that on the upper part there is space enough (ca. 30 cm) to allow the ascent/descent movements of the cover. Regulate the 4 rubber feet to get a perfect horizontal position. Connect the instrument to the power supply socket (230V~ 50/60Hz) using the provided cable.

Important:

- Do not use any extension or adapter and do not modify the provided cable.
- Respect the following working environmental conditions:

Room temperature between 10°C and 40°C.

Relative humidity between 30% and 75%.

Maximum temperature inlet water: 40°C.

3) WATER CONNECTIONS

Inlet running water

Connect the provided pipe (Ø 10 mm) between the connector placed on the left side of the instrument and a water tap.

Outlet running water


Connect the provided pipe (Ø 18 mm) between the connector placed on the left side of the instrument and the drain of the laboratory.

Important: As the draining happens by fall, it is necessary that the pipe is positioned as linear as possible and downhill (the water outlet connector of the stainer must be situated over the water outlet connector of the laboratory).


Connection between the two running water dishes and the instrument

Connect the two provided pipes (Ø 10 mm) between the connectors of the running water dishes and the connectors situated on the right side of the instrument. Choose the position of the dishes according to the main programme. Make sure that the pipes are long enough to rotate the dishes holder disk.

Left side of the instrument.



Right side of the instrument.




Adjustment of the pressure reducer valve

The pressure reducer valve regulates the water flow when the pressure of the plant is irregular or not adjustable.

Instructions for the flow adjustment:

- 1) In main menu touch key MANUAL.
- 2) Enable the wash dishes touching keys OFF by WATER 1 and WATER 2.
- 3) Open abundantly the water.
- 4) Pull the knob till you hear a click and rotate it till you get the right quantity of water (clockwise to increase the flow, anticlockwise to reduce it).
- 5) Finished the adjustment, block the valve by pressing the knob till you hear a click.
- 6) Disable the wash dishes touching keys ON by WATER 1 and WATER 2.
- 7) In order to return to menu, touch exit key positioned low on the right.

Enable the manual functions



Enable the water flow in dish 1 (water 1) and dish 2 (water 2)

4) CONTROL PANEL

Touch screen control panel, intuitive, impermeable IP55, bilingual (Italian and English) for the programming/modification/visualization of the parameters of work and the control of the process of coloration.



Initial screen.



Main menu.

FUNCTIONS	DESCRIPTION
MAINTENANCE	<ul style="list-style-type: none"> • USER MAINTENANCE: Access through password, reserved to Technical Assistance Service (machine parameters). • FILTERS MAINTENANCE: Control/updating of filter's life.
MANUAL	Allows to page to manual functions of machine.
MAP	Allows to page to dishes map and relative contents.
ALARMS	Allows to page to page of alarms.
PROGRAMMES	Allows to page to programmes management.
EXIT	Allows to exit from menu.

5) INSTRUMENT'S START

Switch on the instrument by pressing the green main switch placed on the right side of the instrument.
On display appears the initial screen and a sound signal advises that the machine is working. Effect the zero setting of the machine.

Touch the screen by one of the two small flags (Italian/English) to page to main menu of the functions.

Key ALARMS lightens intermittently and you can hear a sound signal.

Touch key ALARMS to page to alarms page (photo 1).

Touch key RESET to zero the machine. The basket holder wheel will position on starting position (0) (photo 2 and 3).

Touch exit key, positioned under key RESET, to exit from page alarms.

Touch key ESC low on the right to return to main menu.



PHOTO 1




PHOTO 2



PHOTO 3

5.1 MANUAL OPERATIONS

In manual function it is possible to manage manually all machine movements.



1 = Direction arrows for manual management of the vertical movements of the baskets.

2 = Visualizes the basket position on the dishes wheel.

3 = Activates manually the function of basket immersion with relative oscillation inside a dish. Touching key OFF/ON is activated the function. To deactivate the immersion function, touch key OFF/ON again.

4 = Manual function for the basket moving of one position on the dishes wheel. Touching key OFF/ON once, the basket will move of only one position. Keeping pressed key OFF/ON, the basket will continue moving from dish to dish.

5 = Touching key OFF/ON is activated the water flow in dish WATER 1.

6 = Touching key OFF/ON is activated water flow in dish WATER 2.

7 = Function that activates manually the basket dripping. Touch key OFF/ON to activate/deactivate the function.

8 = Touching key OFF/ON in any position on the dishes wheel, the basket will return to position 0.

9 = Touching this key, you exit from present screen.

ATTENTION: It isn't possible to exit from manual management of movements if the machine isn't in position of 0.

6) SETUP OF DATE AND HOUR

It is possible to set up date and hour both from the main screen and from the functions menu touching directly key of date and hour (photo 4).

A message appears for confirmation of the modification, confirm touching key YES or exit touching key NO.

If you continue touching key YES, appears the screen with set date and hour.

Set up the desired date and hour and confirm all the set data touching key "CONFIRMS DATE AND HOUR".

Press exit key low on the right to return to menu.



PHOTO 4



PHOTO 5



PHOTO 6

7) MAP OF THE DISHES AND RELATIVE CONTENTS

Touch key MAP in main menu of the functions.

It is possible to insert for each dish: name, date of load and minutes of use of each reagent up to 23 reagents.

The first column on the left shows the position of the dish.

The second column shows the name of the reagent.

The third column shows date and hour of insertion.

The fourth column shows the duration in minutes of the reagent.

In fifth column the square lightens intermittently when the minutes of use of each single reagent are finished.

Touching the corresponding square, it is possible to name the reagent and insert date and duration.



8) PARAMETERS SETUP AND PROGRAMMING

8.0 PARAMETERS SETUP AND COLORATION STEPS

Page to menu of parameters setup and programming, touching key PROGRAMMES in main menu.

Appears the screen of programme.



Proceed with the steps programming touching key STEPS.

Appears the screen that shows 30 squares, one for each step to program.




Starting from the first key P1 (first step of programme) up on the left, set up the desired parameters for each step of programme:

- Position n° = number of desired dish you can visualize in MAP.
- Immersion T. = time of immersion of the basket in one dish, measured in minutes and seconds.
- Dripping T. = shows the dripping time of a basket after the immersion phase outside the dish, measured in minutes and seconds.
- E.V. water 1 = activates (ON) or deactivates (OFF) the flow of running water in dish N. 1 for the whole duration of immersion of the basket.
- E.V. water 2 = activates (ON) or deactivates (OFF) the flow of running water in dish N. 2 for the whole duration of immersion of the basket.



Set up the parameters in each step, touching on the screen the value to set up.


Appears the screen where it is possible to set up or change the value, type in the value and press ENTER.




The value will be visualized in the step page.
Touch the exit key low on the right to pass to the programming of the following step.
Proceed with the operation above described for each step.

8.1 PROGRAMME SAVING

Once set up the parameters in the relative steps, it is necessary to name and save the programme.
Exit from the screen of steps setup touching the exit key low on the right.
In order to save the steps and relative set up parameters, touch key PROGRAMMES up on the right.
Appears the screen where it is possible to manage the realized programmes.



The steps and relative inserted parameters, are visualized in the square MACHINE STEPS.
Transfer the visualized programme from square MACHINE STEPS to programmes list PROGRAMMES STEPS touching the key LOAD FROM MACHINE.



Touch the number 0. to open the column of programmes.



Select the numbered square where you want to save the realized programme and touch the key SAVE PROGRAMME. In order to name or name again an existing programme, touch key PROGRAMME NAME. Touch the programme number to name or name again in the relative square, write the desired name and press ENTER.



Press again key PROGRAMME NAME in order to disable this function.
The realized and saved programme is so saved and can be retrieved at any time.

8.2 RETRIEVAL OF A PROGRAMME FROM MEMORY

Touch key PROGRAMMES in main menu.

Touch key PROGRAMMES up on the right in the following screen.

Touch key 0. to visualize the column of programmes and use the direction arrows to look for the programme to execute.

Select the programme to execute touching number/name of the programme.

In the first line up appear number and name of the selected programme.

Touch key LOAD PROGRAMME.

In square STEPS PROGRAMMES appear the steps of the retrieved programme.


To execute the retrieved programme, touch key TRANSFER IN MACHINE.

At this point the machine is ready to execute the programme.

Press ESC to exit.

Press START CYCLE to start the programme.

8.3 PROGRAMMES IN EXECUTION





During the execution of the programme, the movements of the machine are visualized in real time.

If in the programme is expected the use of the dishes with running water, the machine activates automatically, before the start of the programme, the water flow for 25 seconds, time necessary for the total replacement of water in the dishes. If the operator wants however to exclude the use of water, also if expected in the cycle, he can do it deactivating the function WATER 1 and WATER 2 touching the corresponding key ON/OFF. The running water flow in the dishes will last as long as the time of immersion of the basket.


When the basket will move to the following station after the wash, the machine will activate again the water flow for 25 seconds so that there will be clean water to the following passage of the basket.

It is, moreover, possible to change the rotation and oscillation speed of the basket touching key SPEED and dragging the cursor on the dial in the desired position.



If necessary, it is possible to interrupt temporarily the programme touching key PAUSE. The message CYCLE IN PAUSE AFTER DRIPPING will appear, the machine stops after the dripping phase of the basket.

To restart the programme, touch again key PAUSE.



If you want to interrupt definitively the coloration cycle, press the emergency button on the frontal panel and, subsequently, effect a reset of the drives (see paragraph 5) START OF THE INSTRUMENT).

8.4 CYCLE END

When the basket will effect the last immersion expected from the programme, the machine will stop the basket in phase of immersion uttering a sound and visual signal.

Touching key START CYCLE the machine will move the basket to load/discharge station 0.

In this way it will be possible to load and unload the baskets always from the machine frontal side.

8.5 PREDEFINED PROGRAMMES

Three predefined coloration protocols from number 10 to number 12 are present in the machine memory in the column of programmes.

In succession, you find schedules relevant to programmes inserted where are shown the programme steps and relative schedule of the reagents map to allow an optimum use of the programmes.

PROGRAMME 10

For the preparation of Autocolor for histology (Hematoxylin Eosin).

Step	Reagent	Code	Programme 11 histology (Hematoxylin Eosin)
1	Bio Clear	06-1782Q	5 minutes
2	Bio Clear		5 minutes
3	Dehyol absolute	06-10077Q	2 minutes
4	Dehyol absolute		2 minutes
5	Dehyol 95	06-10070Q	2 minutes
6	Distilled water		4 minutes
7	Mayer Ematoxylin	05-06002/L	5 minutes
8	Running water		5 minutes
9	Aqueous eosin	05-10002/L	5 minutes
10	Running water		5 minutes
11	Dehyol 95	06-10070Q	30 seconds
12	Dehyol absolute	06-10077Q	1 minute
13	Dehyol absolute		1 minute
14	Bio Clear	06-1782Q	1 minute
15	Bio Clear		3 minutes

PROGRAMME 11

For the preparation of Autocolor for cytology (Papanicolau).

Step	Reagent	Code	Programme 12 cytology (Papanicolau)
1	Dehyol 95	06-10070Q	1 minute
2	Dehyol 95		1 minute
3	Distilled water		2 minutes
4	Harris hematoxylin	05-12011/L	1 minute
5	Running water		5 minutes
6	Dehyol 95	06-10070Q	15 seconds
7	OG6	05-12013/L	2 minutes
8	Dehyol 95	06-10070Q	15 seconds
9	EA50	05-12019/L	5 minutes
10	Dehyol 95	06-10070Q	15 seconds
11	Dehyol absolute	06-10077Q	30 seconds
12	Dehyol absolute		1 minute
13	Bio Clear	06-1782Q	1 minute
14	Bio Clear		2 minutes
15	Bio Clear		2 minutes

In the following schedule you find the preparation of the dishes on Autocolor indistinctly for the execution of programmes 10 and 11.

Position of the reagents to insert in machine

Dish	Reagent	Code
1	Bio Clear	06-1782Q
2	Bio Clear	
3	Dehyol absolute	06-10077Q
4	Dehyol absolute	
5	Dehyol 95	06-10070Q
6	Distilled water	
7	Mayer Ematoxylin	05-06002/L
8	Running water	
9	Aqueous eosin	05-10002/L
10	Running water	
11	Dehyol 95	06-10070Q
12	Dehyol absolute	06-10077Q
13	Dehyol absolute	
14	Bio Clear	06-1782Q
15	Bio Clear	
16	Distilled water	
17	Harris hematoxylin	05-12011/L
18	Dehyol 95	06-10070Q
19	OG6	05-12013/L
20	Dehyol 95	06-10070Q
21	EA50	05-12019/L

PROGRAMME 12

For the preparation of Autocolor for cytology (Papanicolaou) with two baskets at the same time.

	Reagent	Code	Programme 12 cytology (Papanicolaou)
1	Dehyol 95	06-10070Q	2 minutes (*)
2	Distilled water		2 minutes
3	Harris hematoxylin	05-12011/L	1 minute
4	Running water		5 minutes
5	Dehyol 95	06-10070Q	15 seconds
6	OG6	05-12013/L	2 minutes
7	Dehyol 95	06-10070Q	15 seconds
8	EA50	05-12019/L	5 minutes
9	Dehyol 95	06-10070Q	15 seconds
10	Dehyol absolute	06-10077Q	30 seconds
11	Dehyol absolute	06-10077Q	1 minute
12	Bio Clear	06-1782Q	3 minutes
24	Bio Clear	06-1782Q	0 minutes

(*) and (**) baskets start

In the following schedule you find the preparation of Autocolor for the execution of programme of cytology (Papanicolau) with two baskets at the same time.

Position of the reagents to insert in machine

Dish	Reagent	Code
0	Bio Clear	06-1782Q
1	Bio Clear	
2	Dehyol 95	06-10070Q
3	Dehyol 95	
4	Distilled water	
5	Distilled water	
6	Harris hematoxylin	05-12011/L
7	Harris hematoxylin	
8	Running water	
9	Running water	
10	Dehyol 95	06-10070Q
11	Dehyol 95	
12	OG6	05-12013/L
13	OG6	
14	Dehyol 95	06-10070Q
15	Dehyol 95	
16	EA50	05-12019/L
17	EA50	
18	Dehyol 95	06-10070Q
19	Dehyol 95	
20	Dehyol absolute	06-10077Q
21	Dehyol absolute	
22	Dehyol absolute	
23	Dehyol absolute	

9) ALARMS

In main menu is present the key ALARM that allows the access to the page of alarms.

If during the execution of the programme a malfunction should arise, the machine stops and at the same time enters into action a sound and visual signal in the page of programmes in execution.

Touching the flashing key ALL. you enter in the page of noticed alarms.



Check, and if necessary reset, the signalled fault.

The machine will ask to effect two RESET of the alarms. Touch key RESET to reset the alarms; the machine will move in position of vertical safety.

When the machine will reach the position of vertical safety, touch a second time key RESET.

The machine will move in position of zero, ready for the new execution of the programme.

Touch exit key low on the right.



10) MAINTENANCE

It is possible to page to menu maintenance touching key MAINTENANCE in main menu.



Menu MAINTENANCE is divided into two distinct parts (photo 1):

1) The first, MAINTENANCE USER, is accessible only through password and allows to qualified technicians to check all the machine parameters (photo 2 and 3).



PHOTO 1



PHOTO 2





PHOTO 3

2) The second, MAINTENANCE FILTERS, allows the operator to check and update the filters hourcounter in case of filter's replacement. To page to filters maintenance, touch key MAINTENANCE FILTERS.
It is possible to zero/modify the hourcounter by pressing on the shown value.



11) FILTERS

The instrument is provided with an activated charcoal cartridge filter specific for alcohols and xilol, placed in the back side of the instrument. The panel is provided with a digital hourcounter preset at 600 hours (middle life of a filter). For the visualization of the filters hourcounter see the preceding paragraph. When the hourcounter reaches the value of zero hours, an acoustic and visual alarm enters into action in order to advise the operator that the filter must be replaced.

In order to proceed to the replacement, unscrew the two knobs on the left side, remove the cover, extract the exhausted filter and replace it with a new one.

Enter in menu MAINTENANCE FILTERS (see preceding paragraph) and update the filters hourcounter.

Important:

- Before effecting the replacement, switch off the instrument and unplug the power supply cable from the socket.
- In case of accidental infiltration of liquids or foreign bodies in the filter area, contact immediately the Technical Assistance Service and don't use the instrument.
- Dispose the exhausted filter in compliance with the regulation in force.

12) MAIN SAFETY ELEMENTS


- 1) If necessary, it is possible to stop the machine movements by pressing the red emergency button placed in the frontal side of the instrument. In order to restart the instrument, you must rotate the emergency button clockwise and effect a reset of the drives (see paragraph 9) ALARMS).
- 2) The machine is provided with a sensor that detects the position of the pvc rotating table. If the table goes accidentally out of position, the machine stops all the movements and signals the alarm "table out of position".
- 3) A float sensor is positioned inside the discharge central sink. If there is a malfunction of the electric valves or an obstruction of the discharge pipes that prevent the water from flowing, with consequent raising of the level, the machine interrupts the water flow intervening on the electric valves. In this case the machine doesn't stop but signals the fault.

When the level returns to normal, automatically the electric valves start to work normally again.

13) INSTRUCTIONS FOR THE ROTATION OF THE DISH HOLDER DISK

The swivel dish holder disk allows to accede easily to the dishes positioned on the back.

- 1) In main menu touch key MANUAL.
- 2) Touch the key shown in the photo (arrow up) to bring the basket holder cover at end run in height.



- 3) Pull the release knob and rotate cautiously the disk in the desired position, taking care to the pipes connected to the dishes as shown in the sequence of photos.



PHOTO 1



PHOTO 2




PHOTO 3

- 4) The disk is provided with a concave seat by the release knob. Before restarting the instrument, it is necessary to hook again the disk in the proper seat and place the knob in the block original position (photo 1).

14) REPLACING THE FUSES

If the instrument doesn't start, check that the electric cable is properly connected, there is current and the two fuses under the instrument's socket aren't burnt. If necessary, change them with fuses of identic value.



- Fuses holder socket (N. 2 delayed fuses of 4A - T4AH250V).

Instructions for the replacement: Switch off the appliance and unplug the cable from the socket. Make a light pressure on the little plastic carter covering the fuses under the power supply input (if necessary using a little screwdriver), change them and close the carter.

15) CLEANING AND MAINTENANCE

Maintenance for the user

A proper standard of hygiene and the routine maintenance of the laboratory are necessary to keep the instrument in good working conditions. The operator is responsible for the proper cleaning of the stainer. In particular, it is necessary to control the state of water inlet/outlet pipes, draining central basin and running water dishes, which must be free of mould, encrustations, etc. which might obstruct the water passage.

Important: Do not use acid or abrasive cleaning products in order to avoid ruining the painted surface or other parts (washers, etc.) and wear anti-cut gloves over or under the common gloves to avoid possible accidents. It is possible to use bleach to clean the pipes.

Control the filters wear and tear and effect the replacement if necessary. The stainer is, however, provided with visual and acoustic alarm in case of filter saturation. See paragraph 11) Filters.

Annual maintenance

The annual maintenance must be effected only by qualified and authorized staff. For information contact the Technical Assistance Service.

16) MAIN COMPONENTS

CODE	DESCRIPTION	QUANTITY
16-1610	Reagents basin	22
16-1621	Water inlet pipe	2,5 m
16-1622	Water outlet pipe	2,5 m
16-1623	Water basin/instrument black pipe	3 m
16-1640	Slide holder basket	2
36-A01ES	Emergency button	1
36-POMAUT	Knob	1
36-VOLFIL	Threaded knob	2
36-PLC	Plc	1
36-TSCREEN	Terminal	1
36-ALIM	Feeder	1
36-UPS	Ups	1
36-MOTOREPP	Square fan	2
36-AZION	Drive	2
36-VENT	Fan	3
36-ELET	Electric valve + coil	2
36-SOLEN	Solenoid	2
36-MOLSOL	Solenoid spring	2
36-22-17-7	Running water basin spacer	2
36-138-278	Connector 3/8in 1/4in	1
36-138-313	Connector 3/8in 1/4in	1
36-138-335	Connector 3/8in	2
36-138-408	Connector 3/8in 1/4in 90°	2
36-179-4970	Connector	4
36-312-7164	Floating switch	1
36-3114	Water female terminal	2
36-4002N	PVC big disc	1
36-4002NP	PVC little disc	1
36-4004	Central basin	1
36-4005	Central basin support column	3
36-4006	External cover	1
36-4007	Basket holder cover	1
36-4008	Delrin pin	48
36-4009	Basket holder track	25
36-4010	Little cover	24
36-4012	Water discharge connector	1
36-4012A	Basin discharge connector	1
36-4013	Little basin discharge connector	2
36-4014	Connection block	1
36-4015	Pinion	1
36-4015A	Black pinion	1
36-4016	Gear M 1.5 Z 160	1
36-4017	Support	1
36-4018N	Plate 600x600	1
36-4019	Worm screw 240 MM	1
36-4020	Spit nut	1
36-4022	Sensor support	1
36-4024	Toothed pulley 18XL037	1
36-4025	Toothed belt 110XL037	1
36-4028	Up/down motor support	1

36-4032	Bearing + fifth wheel 745	2
36-4037	Nut M 45X1,5	2
36-4038	Bearing 627	2
36-4039	Centring pin	1
36-4043	Column 12 H 6X200	2
36-4050	Proximity sensor	5
36-4072	Bearing KH 12/28	2
36-4073	Bearing KH 10/26	1
36-4074	Bearing AX15.28	1
36-4087	Anti-vibrating R2L50 50X16X10	4
36-4100	Worked dish	2
36-4105	Microgovernor for water	1
36-4582	Battery	2
36-6012	Pinion support	1
36-6013	Rotating pinion	1
36-6014	Fixed rubber dinghy	1
36-6015	Central mast	1
36-6016	Triangle with stop	1
36-6016A	Bolt for triangle	1
36-6017	Triangle	3
36-6023	Aluminium pulley 72 MXL 02	1
36-6024	Aluminium pulley 36 MXL 025	1
36-6060ST	Stirrup for battery	1
36-9226	Protection grill X RQ16	3
36-40311	Cover fixing pan	1
36-40312	Perforated mast Ø 14 M8	1
36-40313	Threaded mast Ø 12	1
36-40455N	Panel	1
37-AR09E2F6A	6 ampere wire filter	1
37-CFA	Active charcoal filter	1
37-INN1/4	1/4 quick connection	6
37-MANTMB	Steel handle	1
37-129-432	Green waterproof switch	1
37-622	Green bulb	1
37-1007	Microswitch	6
37-1045	Power supply cable	1
37-4050F	Female connector 4/50	2
MU16-1600/T	User manual	1

17) PROBLEM SOLVING

The instrument doesn't start

- 1) Verify the correct insertion of the power supply cable and the presence of current in the electric system.
- 2) Control that the two fuses placed under the instrument's socket aren't burnt.
- 3) If the problem continues, contact the Technical Assistance Service.


The instrument start but water doesn't arrive to dishes


- 1) Verify the correct adjustment of the reducing valve (see paragraph 3) WATER CONNECTIONS).
- 2) Verify that the electric valves are enabled (see paragraph 3) WATER CONNECTIONS).
- 3) If the problem continues, contact the Technical Assistance Service.


Water leak from dishes or from central basin

- 1) Control that the outlet pipe is correctly positioned (the outlet connector of the stainer must be situated over the outlet connector of the laboratory).
- 2) Control that the pressure reducer valve is correctly regulated.
- 3) Verify the seal of inlet/outlet connectors of the running water dishes.
- 4) Control that internal and external outlet pipes aren't blocked. If necessary, effect the cleaning using bleach or similar products.
- 5) Check the correct working of the emergency float.

18) ELECTRIC WIRING DIAGRAM







Milano, 18 giugno 2013

13) DICHIARAZIONE DI CONFORMITA' / DECLARATION OF CONFORMITY KONFORMITÄTSEKTLÄRUNG / DECLARATION DE CONFORMITE

Nome e indirizzo della ditta Name and address of the firm Name und Adresse der Firma Nom et adresse de l'entreprise		BIO OPTICA Milano S.p.A. Via S.Faustino, 58 20134 MILANO C.F./P.IVA 06754140157
<p align="center">Dichiariamo sotto la nostra responsabilità che / We declare under our sole responsibility that Wir erklären in alleiniger Verantwortung, dass / Nous déclarons sous notre propre responsabilité que</p>		
il dispositivo medico-diagnostico in vitro the in vitro diagnostic medical device das Medizinprodukt für die In-vitro_Diagnostik le dispositif médical de diagnostic in vitro		AUTOCOLOR TOUCH
della classe: of class: der Klasse: de la classe:	Altro Other Sonstiges produkti Autre	
<p align="center">soddisfa tutte le disposizioni della direttiva 98/79/CE e successive modifiche ed integrazioni che lo riguardano meets all the provisions of the directive 98/79/EC and following amendment which apply to it allen Anforderungen der Richtlinie 98/79/EG entspricht, die anwendbar sind remplit toutes les exigences de la directive 98/79/CE et modification qui lui sont applicables</p>		
Norme nazionali o armonizzate applicate Applied harmonised standards and National standards Angewandte harmonisierte Normen, nationale Normen Normes harmonisées et normes nationales		– EN 375 – EN 980 – EN ISO 14971 – EN 60601-1-2 – EN 61010-1 – EN 61010-2-101
		Edizione in vigore alla data di emissione del documento/Current ed. at document date/Aktuelle Ausgabe am belegdatum/Édition actuelle à la date du document

BIO-OPTICA MILANO SPA
Legale Rappresentante
Carlo Sbona

